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Instructions for use BIOLAN®

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Eye drops with sodium hyaluronate 0.15 %

DESCRIPTION

BIOLAN® form a sterile, preservativefree, elastoviscous, clear protective film for the cornea. **BIOLAN®** contain Sodium Hyaluronate, the sodium salt of the Hyaluronic acid, which forms, due to its physical properties, an evenly, solid and long-retaining protective shield to the eye's epithelial layer without impairing vision. As a result, **BIOLAN®** cause long-lasting relief from discomforts associated with: • mechanical stress, e.g. during wearing of hard or soft contact lenses, or during ophthalmic diagnostic procedures • environmental stress, e.g. induced by air conditioning, wind, coldness, dryness or air polluted by cigarette smoke for example • visual stress from intensive viewing, e.g. at computer terminals, microscopes or during prolonged operation of motor vehicles. Preservative free **BIOLAN®** are also well-tolerated during long-term use. **BIOLAN®** are suitable for contact lenses wearers. The contact lenses do not have to be removed during the application of the eye drops.

APPLICATION

Detach one single dose container. Make sure that each single dose container is intact. To open the single dose container, twist off the flat end (do not pull). Tilt your head back a little, look upwards and gently pull the lower eyelid away from the eye and instill one drop in the conjunctival sac of the eye 3 - 5 times daily. Close your eyes slowly to allow an even distribution of the liquid on the eye surface. If the discomfort gets worse or does not improve, please consult your ophthalmologist!

COMPOSITION

1 ml contains: sodium hyaluronate 0.15 %, sodium chloride, sodium citrate, citric acid, water for injections.

PACKAGE SIZES

Commercial units: single dose containers, 10 x 0.35 ml / 20 x 0.35 ml / 40 x 0.35 ml / 60 x 0.35 ml

PRECAUTIONS AND WARNINGS

- Do not use the product if single dose container is damaged. • Do not touch the eye with the single dose container.
- Do not use if you are allergic to any of the ingredients. • For single use only; as **BIOLAN®** do not contain preservatives, any solution not used immediately after opening, should be discarded. • If you also use other eye drops / eye ointment, wait about 15 minutes. **BIOLAN®** should always be used last.
- Do not use in infected or injured eyes.
- Temporary blurring for a brief period may be experienced after application; avoid driving, operating machinery and wait until normal vision is restored. • Do not use the eye drops beyond the expiry date which you can find on each single dose container and on the folding box.

STORAGE

Store at 2 °C to 25 °C. Do not freeze. Protect from light by storing the single dose containers in the folding box.

Lot-number and expiry date are embossed on each single dose container.

MANUFACTURER:
Penta Arzneimittel GmbH
Werksstrasse 3
92551 Stulln
Germany

DISTRIBUTION IN SWITZERLAND:

medServe AG
Churerstrasse 35
8808 Pfäffikon SZ
Switzerland

for single use only.

aseptic manufacturing.

For description and mode of application see attached instructions for use.

storage temperature.

Do not use the product if single dose container is damaged.

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