




an-Bandage Lenses



A transparent protection to cover corneal lesions resulting from abrasions and erosions or corneal ulcers

Technical Data:

Central Thickness	0.3 mm
Power	0.0 – 0.25 D
Material	Polyxylon
Water Content	74%
Dk (oxygen)	33

- Exact fitting lenses with a high retention rate
- Immediate pain relief in case of corneal lesion
- Designed specifically for animals
- The small “a” indicates the correct fit (not inside out)

an-vision GmbH

Neuendorfstr 22a
16761 Hennigsdorf
Germany
T: +49 (0)3302-209 5550
info@an-vision.eu

www.an-vision.com

an-vision Inc.

9067 S.1300 W.
Salt Lake City, UT 84088
USA
T: 801-561-5040
info@an-vision.com

Sizes of an-Bandage Lenses

Cats

	Diameter	Base Curve	
C1	16.0 mm	9.2 mm	Cats normal - slightly larger
C2	14.0 mm	8.3 mm	Kitten
C3	15.5 mm	8.5 mm	Young Cats
C4	16.0 mm	8.8 mm	Cats normal - slightly smaller
C5	18.0 mm	9.4 mm	Persian Cats - large eyes

Dogs

	Diameter	Base Curve	
D1	18.0 mm	9.6 mm	Medium sized dogs - somewhat smaller; most dog breeds such as collies, huskies, corgis, some boxers
D2	18.0 mm	9.8 mm	Medium sized dogs - somewhat larger; most dog breeds such labradors, retrievers, boxers
D3	20.0 mm	10.4 mm	Dogs - large eyes; some bulldogs; some boxers
D4	15.5 mm	8.5 mm	Young dogs
D5	22.0 mm	11.4 mm	Dogs - very large eyes
D6	22.0 mm	11.8 mm	Dogs - very large eyes only
D7	20.0 mm	11.2 mm	Dogs - in between sizes
D8	20.0 mm	10.8 mm	Dogs - in between sizes
D9	18.0 mm	10.0 mm	Dogs - in between sizes

Horses

	Diameter	Base Curve	
H1	34.0 mm	18.0 mm	Horses, Normal Eye Size
H2	32.0 mm	17.0 mm	Small Horses, Ponies
H3	26.0 mm	15.0 mm	Young Ponies
H4	30.0 mm	16.5 mm	Young Horses
H5	38.0 mm	19.5mm	Large Horses
H6	36.0 mm	20.0 mm	Horses, Larger Eyes