



## an- Test Strips

### Ophthalmic Diagnostic Test Strips



Measurement of tear production  
Each strip is marked with a measuring scale and impregnated with a blue or white indicator



Staining dead and degenerated cells and mucus.  
Diagnose of KSC and superficial corneal epithelial abnormalities. Each strip is impregnated with 1.5 mg Rose Bengal



To detect corneal and conjunctival defects and to assess nasolacrimal patency.  
Each strip is impregnated with Fluorescein Sodium



Staining of dead or degenerated conjunctival cells  
Early detection of dry eye  
Each strip is impregnated with 1.5 mg Lissamine Green

**an-vision GmbH**  
Neuendorfstr 22a  
16761 Hennigsdorf  
Germany  
T: +49 (0)3302-209 5550  
[info@an-vision.eu](mailto:info@an-vision.eu)

[www.an-vision.com](http://www.an-vision.com)

**an-vision Inc.**  
9067 S.1300 W.  
Salt Lake City, UT 84088  
USA  
T: 801-561-5040  
[info@an-vision.com](mailto:info@an-vision.com)